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# haptosil D Knetsilikon

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## Instructions for use

# haptosil D – kneading silicone

Before using the product, please read these instructions for use carefully.

### 1. Product description

haptosil D is an addition-curing kneading silicone with a Shore A hardness of 90 for the use in dental laboratories.

### 2. Indication range

haptosil D is used for the fabrication of keys, repair models and checkbits.

### 3. Safety recommendations and hazard warnings

No special measures required. After eye contact, rinse eyes under running water for several minutes with lid being opened.

If complaints persist, seek medical attention.

### 4. Processing

Use the enclosed dosing spoons to remove the two components from the jars. Make sure not to mix up dosing spoons and lids when removing material from the jars! haptosil D base and catalyst paste are mixed in the ratio of 1:1 and the paste is mixed until it is completely free from streaks (approx. 30 sec.). Make sure to adhere exactly to the given mixing ratio to achieve the required setting times and final hardness. The jars must be sealed after the use.

It is recommended to wear vinyl gloves when mixing the pastes to avoid contamination. Store at room temperature. Higher temperatures will reduce the processing and setting time and lower temperatures will increase them. Materials containing sulfur or amine may result in retarded setting or even prevent setting (e.g. various milling oils, latex gloves).

### 5. Technical data (times from beginning of mixing, calculated at 23° C)

Mixing time:	30 seconds
Processing time:	60 seconds
Setting time:	7 minutes
Shore A hardness (after 24 hrs):	≥ 90 Shore A

### 6. Disposal

The hardened silicone can be disposed of together with commercial or household waste.

### 7. Additional information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend to read the instructions for use again before using a new package.

Any recommendations concerning the use are based on brendent's experience and tests, and can only be considered to provide approximate values. Our products are subject to continuous further development. Therefore we reserve the right to perform changes in the design and composition.

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