


WORKFLOW DIAGRAM



CUTTING GUIDE

Non-engaging

Platform	Height (mm)	Reference
Narrow	3.5	IPD/KA-IN-35
	5.0	IPD/KA-IN-50
	6.5	IPD/KA-IN-65
Standard	3.5	IPD/KA-IR-35
	5.0	IPD/KA-IR-50
	6.5	IPD/KA-IR-65
Wide	3.5	IPD/KA-IW-35
	5.0	IPD/KA-IW-50
	6.5	IPD/KA-IW-65

Version with agulation window subject to availability.

Engaging



Platform	Height (mm)	Reference
Narrow	3.5	IPD/KR-IN-35
	5.0	IPD/KR-IN-50
	6.5	IPD/KR-IN-65
Standard	3.5	IPD/KR-IR-35
	5.0	IPD/KR-IR-50
	6.5	IPD/KR-IR-65
Wide	3.5	IPD/KR-IW-35
	5.0	IPD/KR-IW-50
	6.5	IPD/KR-IW-65

Version with agulation window subject to availability.

FIXING SCREW

	Screwdriver	Reference
Every cutting guide	1.20 Hex	IPD/KA-TI-00

CUTTING GUIDE KITS

Platform	Heights (mm)		
Narrow	3.5 / 5.0 / 6.5	IPD/KT-IN-00	IPD/KT-IN-01
Standard	3.5 / 5.0 / 6.5	IPD/KT-IR-00	IPD/KT-IR-01
Wide	3.5 / 5.0 / 6.5	IPD/KT-IW-00	IPD/KT-IW-01

All kits include the locking screw IPD/KA-TI-00

SCAN BODY SCREWDRIVER

	Scan Body height (mm)	Reference
	10	IPD/PH-20-18
	15	IPD/PH-20-25



Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides
Nobel Biocare® Branemark System® Klockner® KL™	Ø3.5	Engaging	1.5	IPD/AA-IN-02/3D	IPD/KT-IN-00
		Non-engaging		IPD/AA-IN-03/3D	IPD/KT-IN-01
Nobel Biocare® Branemark System® Biomet® 3i® Externa™ 3i® BTI® Hex. Ext. Sweden & Martina® Outlink® Microdent® Microdent® Universal™ Klockner® KL™	Ø4.1	Engaging	0.5	IPD/AA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/AA-IR-02/3D	
		Non-engaging	2.5	IPD/AA-IR-04/3D	IPD/KT-IR-01
			0.5	IPD/AA-IR-01/3D	
Phibo® TSH®	S3-S4	Non-engaging	1.5	IPD/AA-IR-03/3D	IPD/KT-IR-01
Nobel Biocare® Branemark System®	Ø5.1	Engaging	1.3	IPD/AA-IW-02/3D	IPD/KT-IW-00
		Non-engaging		IPD/AA-IW-03/3D	-
Nobel Biocare® Multii-unit System	Ø4.8	Engaging	0.5	IPD/AB-IR-00/3D	-
			0.5	IPD/AB-IR-02/3D	-
		Non-engaging	0.2	IPD/AB-IR-01/3D	-
			0.2	IPD/AB-IR-03/3D	-
Nobel Biocare® Replace® Select	Ø3.5	Engaging	0.5	IPD/AC-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/AC-IN-02/3D	
			2.5	IPD/AC-IN-04/3D	
	Non-engaging	0.5	IPD/AC-IN-01/3D	IPD/KT-IN-01	
		1.5	IPD/AC-IN-03/3D		
		2.5	IPD/AC-IN-05/3D		
Ø4.3	Engaging	0.5	IPD/AC-IR-00/3D	IPD/KT-IR-00	
		1.5	IPD/AC-IR-02/3D		
	Non-engaging	0.5	IPD/AC-IR-01/3D	IPD/KT-IR-01	
		1.5	IPD/AC-IR-03/3D		
Ø5.0	Engaging	0.5	IPD/AC-IW-00/3D	IPD/KT-IW-00	
	Non-engaging	1.5	IPD/AC-IW-02/3D	IPD/KT-IW-01	
Ø5.0	Non-engaging	0.5	IPD/AC-IW-01/3D	IPD/KT-IW-01	
		1.5	IPD/AC-IW-03/3D		

Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides
Nobel Biocare® Nobel Active®	Ø3.0	Engaging	0.8	IPD/AD-IT-00/3D	IPD/KT-IN-00
			1.5	IPD/AD-IT-02/3D	
			2.5	IPD/AD-IT-04/3D	
		Non-engaging	0.8	IPD/AD-IT-01/3D	IPD/KT-IN-01
			1.5	IPD/AD-IT-03/3D	
	Ø3.5	Engaging	2.5	IPD/AD-IT-05/3D	IPD/KT-IN-00
			0.5	IPD/AD-IN-00/3D	
			1.5	IPD/AD-IN-02/3D	
		Non-engaging	2.5	IPD/AD-IN-04/3D	IPD/KT-IN-01
			3.5	IPD/AD-IN-06/3D	
Biomet® 3i® Externa™ 3i®	Ø4.3	Engaging	0.5	IPD/AD-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/AD-IN-03/3D	
			2.5	IPD/AD-IN-05/3D	
		Non-engaging	3.5	IPD/AD-IN-07/3D	IPD/KT-IR-00
			0.5	IPD/AD-IR-00/3D	
	Ø3.4	Engaging	1.5	IPD/AD-IR-02/3D	IPD/KT-IR-01
			2.5	IPD/AD-IR-04/3D	
			0.5	IPD/AD-IR-01/3D	
		Non-engaging	1.5	IPD/AD-IR-03/3D	IPD/KT-IR-00
			2.5	IPD/AD-IR-05/3D	
Microdent® Microdent® Universal™	Ø5.0	Engaging	0.5	IPD/BA-IN-02/3D	IPD/KT-IN-00
			1.5	IPD/BA-IN-03/3D	IPD/KT-IN-01
	Non-engaging	0.5	IPD/BA-IW-00/3D	IPD/KT-IW-00	
		1.5	IPD/BA-IW-02/3D	IPD/KT-IW-01	
Biomet® 3i® Certain®	Ø5.1	Engaging	0.5	IPD/BA-IW-01/3D	IPD/KT-IW-01
			1.5	IPD/BA-IW-03/3D	
			0.5	IPD/BB-IN-00/3D	
		Non-engaging	1.5	IPD/BB-IN-02/3D	IPD/KT-IN-00
			2.5	IPD/BB-IN-04/3D	
	Ø3.4	Engaging	0.5	IPD/BB-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/BB-IN-03/3D	
			2.5	IPD/BB-IN-05/3D	
		Non-engaging	0.5	IPD/BB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/BB-IR-02/3D	IPD/KT-IR-01
Ø4.1	Engaging	0.5	IPD/BB-IR-01/3D	IPD/KT-IR-01	
		1.5	IPD/BB-IR-03/3D		
		0.5	IPD/BB-IW-00/3D		
	Non-engaging	1.5	IPD/BB-IW-02/3D	IPD/KT-IW-00	
		0.5	IPD/BB-IW-01/3D	IPD/KT-IW-01	
Ø5.0	Non-engaging	1.5	IPD/BB-IW-03/3D	IPD/KT-IW-01	

Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides	
Klockner® SK2-NK2	Ø4.2	Engaging	2.25	IPD/CA-IR-04/3D	IPD/KT-IR-00	
		Non-engaging		IPD/CA-IR-05/3D	IPD/KT-IR-01	
Klockner® Essential Cone®	Ø4.5	Engaging	0.5	IPD/CB-IR-00/3D	IPD/KT-IR-00	
		Non-engaging		IPD/CB-IR-01/3D	IPD/KT-IR-01	
		Engaging		0.7	IPD/CD-IN-00/3D	IPD/KT-IN-00
				1.5	IPD/CD-IN-02/3D	
Ø3.5		Engaging	2.5	IPD/CD-IN-04/3D	IPD/KT-IN-01	
			0.7	IPD/CD-IN-01/3D		
		Non-engaging	1.5	IPD/CD-IN-03/3D	IPD/KT-IN-01	
			2.5	IPD/CD-IN-05/3D		
Klockner® Vega®	Ø4.1/Ø4.5	Engaging	0.7	IPD/CD-IR-00/3D	IPD/KT-IR-00	
			1.5	IPD/CD-IR-02/3D		
			2.5	IPD/CD-IR-04/3D		
		Non-engaging	0.7	IPD/CD-IR-01/3D	IPD/KT-IR-01	
			1.5	IPD/CD-IR-03/3D		
			2.5	IPD/CD-IR-05/3D		
Straumann® Tissue Level	Ø4.8	Engaging	0.4	IPD/DA-IR-00/3D	IPD/KT-IR-00	
		Non-engaging		IPD/DA-IR-01/3D	IPD/KT-IR-01	
	Ø6.5	Engaging	1.25	IPD/DA-IW-02/3D	IPD/KT-IW-00	
		Non-engaging		IPD/DA-IW-03/3D	IPD/KT-IW-01	
	Ø3.3 NC		Engaging	0.5	IPD/DB-IN-00/3D	IPD/KT-IN-00
				1.5	IPD/DB-IN-02/3D	
2.5				IPD/DB-IN-04/3D		
Non-engaging			3.5	IPD/DB-IN-06/3D	IPD/KT-IN-01	
			0.5	IPD/DB-IN-01/3D		
			1.5	IPD/DB-IN-03/3D		
Straumann® Bone Level		Engaging	2.5	IPD/DB-IN-05/3D	IPD/KT-IR-00	
			3.5	IPD/DB-IN-07/3D		
			0.5	IPD/DB-IR-00/3D		
		Non-engaging	1.5	IPD/DB-IR-02/3D	IPD/KT-IR-00	
			2.5	IPD/DB-IR-04/3D		
			3.5	IPD/DB-IR-06/3D		
Ø4.1 RC		Engaging	0.5	IPD/DB-IR-01/3D	IPD/KT-IR-01	
			1.5	IPD/DB-IR-03/3D		
			2.5	IPD/DB-IR-05/3D		
		Non-engaging	3.5	IPD/DB-IR-07/3D	IPD/KT-IR-01	
			1	IPD/DC-IR-00/3D		
			2	IPD/DC-IR-02/3D		
Straumann® BLX®	RB/WB	Engaging	3	IPD/DC-IR-04/3D	IPD/KT-IR-00	
			1	IPD/DC-IR-01/3D		
			2	IPD/DC-IR-03/3D		
		Non-engaging	3	IPD/DC-IR-05/3D	IPD/KT-IR-01	
			1	IPD/DC-IR-01/3D		

Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides
Astra® Osseospeed®	Ø3.0	Engaging	0.8	IPD/EA-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/EA-IN-02/3D	
			2.5	IPD/EA-IN-04/3D	
		Non-engaging	0.8	IPD/EA-IN-01/3D	IPD/KT-IN-01
			1.5	IPD/EA-IN-03/3D	
			2.5	IPD/EA-IN-05/3D	
	Ø3.5-4.0	Engaging	0.5	IPD/EA-IR-00/3D	IPD/KT-IN-00
			1.5	IPD/EA-IR-02/3D	
			2.5	IPD/EA-IR-04/3D	
		Non-engaging	3.5	IPD/EA-IR-06/3D	IPD/KT-IN-01
			0.5	IPD/EA-IR-01/3D	
			1.5	IPD/EA-IR-03/3D	
Ø4.5-5.0	Engaging	2.5	IPD/EA-IR-05/3D	IPD/KT-IR-00	
		3.5	IPD/EA-IR-07/3D		
		0.5	IPD/EA-IW-00/3D		
	Non-engaging	1.5	IPD/EA-IW-02/3D	IPD/KT-IR-01	
		0.5	IPD/EA-IW-01/3D		
		1.5	IPD/EA-IW-03/3D		
Astra® Evolution®	Ø3.6	Engaging	0.5	IPD/EB-IN-00/3D	IPD/KT-IN-00
			1.5	IPD/EB-IN-02/3D	
			2.5	IPD/EB-IN-04/3D	
	Non-engaging	0.5	IPD/EB-IN-01/3D	IPD/KT-IN-01	
		1.5	IPD/EB-IN-03/3D		
		2.5	IPD/EB-IN-05/3D		
Ø4.2	Engaging	0.5	IPD/EB-IR-00/3D	IPD/KT-IR-00	
		1.5	IPD/EB-IR-02/3D		
		2.5	IPD/EB-IR-04/3D		
	Non-engaging	0.5	IPD/EB-IR-01/3D	IPD/KT-IR-01	
		1.5	IPD/EB-IR-03/3D		
		2.5	IPD/EB-IR-05/3D		
Ø3.5	Engaging	0.5	IPD/FA-IN-00/3D	IPD/KT-IN-00	
		1.5	IPD/FA-IN-02/3D		
		2.5	IPD/FA-IN-04/3D		
	Non-engaging	0.5	IPD/FA-IN-01/3D	IPD/KT-IN-01	
		1.5	IPD/FA-IN-03/3D		
		2.5	IPD/FA-IN-05/3D		
Zimmer® Screw Vent®	Ø4.5	Engaging	0.5	IPD/FA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/FA-IR-02/3D	
			2.5	IPD/FA-IR-04/3D	
	Non-engaging	0.5	IPD/FA-IR-01/3D	IPD/KT-IR-01	
		1.5	IPD/FA-IR-03/3D		
		2.5	IPD/FA-IR-05/3D		
Ø5.7		Engaging	0.5	IPD/FA-IW-00/3D	IPD/KT-IW-00
		Non-engaging	0.5	IPD/FA-IW-01/3D	IPD/KT-IW-01

Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides
Zimmer® Eztetic®	Ø3.1	Engaging	1	IPD/FC-IT-00/3D	IPD/KT-IN-00
			2	IPD/FC-IT-02/3D	
			3	IPD/FC-IT-04/3D	
		Non-engaging	1	IPD/FC-IT-01/3D	
			2	IPD/FC-IT-03/3D	
			3	IPD/FC-IT-05/3D	
BTI® Tiny BTI®	Ø3.5	Engaging	1.5	IPD/GA-IN-02/3D	IPD/KT-IN-00
		Non-engaging		IPD/GA-IN-03/3D	IPD/KT-IN-01
BTI® Interna Universal®	Ø4.1	Engaging	0.5	IPD/GB-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/GB-IR-02/3D	
			2.5	IPD/GB-IR-04/3D	
		Non-engaging	0.5	IPD/GB-IR-01/3D	
			1.5	IPD/GB-IR-03/3D	
BTI® Multi-im®	Ø4.1	Non-engaging	0.3	IPD/GC-IR-01/3D	-
		Engaging	1.5	IPD/IA-IN-02/3D	IPD/KT-IN-00
Dentsply® Xive® Friadent®	Ø3.4		Non-engaging	1.5	IPD/IA-IN-03/3D
		Engaging	0.5	IPD/IA-IR-00/3D	IPD/KT-IR-00
	1.5		IPD/IA-IR-02/3D		
	Non-engaging		0.5	IPD/IA-IR-01/3D	
			1.5	IPD/IA-IR-03/3D	
	Dentsply® Ankylos®	C/X	Engaging	1.5	IPD/IB-IR-02/3D
2.5				IPD/IB-IR-04/3D	
1.5				IPD/IB-IR-03/3D	
Non-engaging			1.5	IPD/IB-IR-05/3D	
			2.5	IPD/IB-IR-05/3D	
Camlog® Camlog®	Ø3.3	Engaging	1.5	IPD/JA-IN-02/3D	IPD/KT-IN-00
		Non-engaging	1.5	IPD/JA-IN-03/3D	IPD/KT-IN-01
	Ø3.8	Engaging	0.5	IPD/JA-IR-00/3D	IPD/KT-IR-00
			1.5	IPD/JA-IR-02/3D	
		Non-engaging	0.5	IPD/JA-IR-01/3D	
	Ø4.3	Engaging	1.5	IPD/JA-IR-03/3D	IPD/KT-IR-01
			0.5	IPD/JA-IW-00/3D	
		Non-engaging	1.5	IPD/JA-IW-02/3D	
			0.5	IPD/JA-IW-01/3D	
		1.5	IPD/JA-IW-03/3D	IPD/KT-IR-01	

Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides		
Sweden & Martina® Outlink®	Ø3.5	Engaging	1.5	IPD/MA-IN-02/3D	IPD/KT-IN-00		
			2.5	IPD/MA-IN-04/3D			
		Non-engaging	1.5	IPD/MA-IN-03/3D			
	Ø5.0	Engaging	2.5	IPD/MA-IN-05/3D		IPD/KT-IN-01	
			0.5	IPD/BA-IW-00/3D			
		Non-engaging	1.5	IPD/BA-IW-02/3D			
Sweden & Martina® Premium™ Kohno®	Ø3.3	Engaging	0.5	IPD/BA-IW-01/3D	IPD/KT-IW-01		
			0.5	IPD/BA-IW-03/3D			
		Non-engaging	1.5	IPD/BA-IW-03/3D			
	Ø3.8	Engaging	0.5	IPD/MB-IT-00/3D	IPD/KT-IN-00		
			1.5	IPD/MB-IT-02/3D			
			0.5	IPD/MB-IT-01/3D			
		Non-engaging	1.5	IPD/MB-IT-03/3D			
			0.7	IPD/MB-IN-00/3D			
			1.5	IPD/MB-IN-02/3D			
Ø4.25	Engaging	0.7	IPD/MB-IN-01/3D	IPD/KT-IN-01			
		1.5	IPD/MB-IN-03/3D				
	Non-engaging	0.7	IPD/MB-IR-00/3D				
Phibo® TSH®	S2	Engaging	1.5	IPD/MB-IR-02/3D	IPD/KT-IR-00		
		Non-engaging	1.5	IPD/MB-IR-03/3D			
	S5	Engaging	0.7	IPD/MB-IR-01/3D			
		Non-engaging	1.5	IPD/MB-IR-03/3D			
	Osstem® TSIII®	Ø3.5 Mini	Engaging	1.5		IPD/NA-IN-02/3D	IPD/KT-IN-00
				1.5		IPD/NA-IN-03/3D	
Non-engaging			1.3	IPD/NA-IW-02/3D			
Ø4.0 Regular		Engaging	1.3	IPD/NA-IW-03/3D	IPD/KT-IW-01		
			1	IPD/OB-IN-00/3D			
			2	IPD/OB-IN-02/3D			
		Non-engaging	3	IPD/OB-IN-04/3D			
			1	IPD/OB-IN-01/3D			
			2	IPD/OB-IN-03/3D			
Biotech® Dental Kontakt®	RP	Engaging	3	IPD/OB-IN-05/3D	IPD/KT-IR-01		
			1	IPD/OB-IR-00/3D			
		Non-engaging	2	IPD/OB-IR-02/3D			
			3	IPD/OB-IR-04/3D			
			1	IPD/OB-IR-01/3D			
Biotech® Dental Kontakt®	RP	Engaging	2	IPD/OB-IR-03/3D	IPD/KT-IR-01		
			3	IPD/OB-IR-05/3D			
		Non-engaging	1	IPD/QB-IR-00/3D			
			2	IPD/QB-IR-02/3D			
			3	IPD/QB-IR-04/3D			
Biotech® Dental Kontakt®	RP	Engaging	1	IPD/QB-IR-01/3D	IPD/KT-IR-01		
			2	IPD/QB-IR-03/3D			
		Non-engaging	3	IPD/QB-IR-05/3D			
			1	IPD/QB-IR-01/3D			
			2	IPD/QB-IR-03/3D			

Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides	
Neodent® Grand Morse™ GM	GM	Engaging	1	IPD/RB-IR-00/3D	IPD/KT-IR-00	
			2	IPD/RB-IR-02/3D		
			3	IPD/RB-IR-04/3D		
		Non-engaging	1	IPD/RB-IR-01/3D		IPD/KT-IR-01
			2	IPD/RB-IR-03/3D		
			3	IPD/RB-IR-05/3D		
	Ø3.25/3.75	Engaging	0.5	IPD/SB-IN-00/3D	IPD/KT-IN-00	
			1.5	IPD/SB-IN-02/3D		
			2.5	IPD/SB-IN-04/3D		
		Non-engaging	0.5	IPD/SB-IN-01/3D		IPD/KT-IN-01
1.5			IPD/SB-IN-03/3D			
2.5			IPD/SB-IN-05/3D			
Bego® Semados®	Ø4.1	Engaging	0.7	IPD/SB-IR-00/3D	IPD/KT-IR-00	
			1.7	IPD/SB-IR-02/3D		
			2.7	IPD/SB-IR-04/3D		
		Non-engaging	0.7	IPD/SB-IR-01/3D		IPD/KT-IR-01
			1.7	IPD/SB-IR-03/3D		
			2.7	IPD/SB-IR-05/3D		
	Ø4.5	Engaging	0.5	IPD/SB-IW-00/3D	IPD/KT-IR-00	
			1.5	IPD/SB-IW-02/3D		
		Non-engaging	2.5	IPD/SB-IW-04/3D		IPD/KT-IR-01
			0.5	IPD/SB-IW-01/3D		
MIS® Seven®	Narrow Ø3.3	Engaging	1.5	IPD/TA-IT-00/3D	IPD/KT-IN-00	
			2.5	IPD/TA-IT-02/3D		
			2.5	IPD/TA-IT-04/3D		
		Non-engaging	0.5	IPD/TA-IT-01/3D		IPD/KT-IN-01
			1.5	IPD/TA-IT-03/3D		
			2.5	IPD/TA-IT-05/3D		
	Standard	Engaging	0.5	IPD/FA-IN-00/3D	IPD/KT-IN-00	
			1.5	IPD/FA-IN-02/3D		
		Non-engaging	2.5	IPD/FA-IN-04/3D		IPD/KT-IN-01
			0.5	IPD/FA-IN-01/3D		
Wide	Engaging	1.5	IPD/FA-IN-03/3D	IPD/KT-IR-00		
		2.5	IPD/FA-IN-05/3D			
		0.5	IPD/FA-IR-00/3D			
	Non-engaging	1.5	IPD/FA-IR-02/3D		IPD/KT-IR-01	
		2.5	IPD/FA-IR-04/3D			
		0.5	IPD/FA-IR-01/3D			
Wide	Engaging	1.5	IPD/FA-IR-03/3D	IPD/KT-IR-01		
		2.5	IPD/FA-IR-05/3D			
	Non-engaging	1.5	IPD/FA-IR-03/3D		IPD/KT-IR-01	
		2.5	IPD/FA-IR-05/3D			

Compatibility	Platform	Connection	Gingival height (mm)	Reference	Cutting guides	
MIS® C1/V3®	Standard	Engaging	0.7	IPD/TB-IR-00/3D	IPD/KT-IR-00	
			1.5	IPD/TB-IR-02/3D		
			2.5	IPD/TB-IR-04/3D		
		Non-engaging	3.5	IPD/TB-IR-06/3D		IPD/KT-IR-01
			0.5	IPD/TB-IR-01/3D		
			1.5	IPD/TB-IR-03/3D		
	NP	Engaging	2.5	IPD/TB-IR-05/3D	IPD/KT-IR-01	
			3.5	IPD/TB-IR-07/3D		
		Non-engaging	1	IPD/UB-IN-00/3D		IPD/KT-IN-00
			2	IPD/UB-IN-02/3D		
DIO® UFII	RP	Engaging	3	IPD/UB-IN-04/3D	IPD/KT-IN-01	
			1	IPD/UB-IN-01/3D		
			2	IPD/UB-IN-03/3D		
		Non-engaging	3	IPD/UB-IN-05/3D		IPD/KT-IR-00
			1	IPD/UB-IR-00/3D		
			2	IPD/UB-IR-02/3D		
	RP	Engaging	3	IPD/UB-IR-04/3D	IPD/KT-IR-01	
			1	IPD/UB-IR-01/3D		
		Non-engaging	2	IPD/UB-IR-03/3D		IPD/KT-IR-01
			3	IPD/UB-IR-05/3D		
Megagen® AnyOne®	RP	Engaging	1	IPD/WA-IR-00/3D	IPD/KT-IR-00	
			2	IPD/WA-IR-02/3D		
			3	IPD/WA-IR-04/3D		
		Non-engaging	1	IPD/WA-IR-01/3D		IPD/KT-IR-01
			2	IPD/WA-IR-03/3D		
			3	IPD/WA-IR-05/3D		
	RP	Engaging	1	IPD/WB-IR-00/3D	IPD/KT-IR-00	
			2	IPD/WB-IR-02/3D		
		Non-engaging	3	IPD/WB-IR-04/3D		IPD/KT-IR-01
			1	IPD/WB-IR-01/3D		
Megagen® AnyRidge®	RP	Engaging	2	IPD/WB-IR-03/3D	IPD/KT-IR-01	
			3	IPD/WB-IR-05/3D		
			1	IPD/WB-IR-00/3D		
		Non-engaging	2	IPD/WB-IR-02/3D		IPD/KT-IR-00
			3	IPD/WB-IR-04/3D		
			1	IPD/WB-IR-01/3D		
	RP	Engaging	2	IPD/WB-IR-03/3D	IPD/KT-IR-01	
			3	IPD/WB-IR-05/3D		
		Non-engaging	1	IPD/XA-IR-00/3D		IPD/KT-IR-00
			2	IPD/XA-IR-02/3D		
Dentium® Implantium® / Superline™	RP	Engaging	3	IPD/XA-IR-04/3D	IPD/KT-IR-01	
			1	IPD/XA-IR-01/3D		
			2	IPD/XA-IR-03/3D		
		Non-engaging	3	IPD/XA-IR-05/3D		IPD/KT-IR-01
			1	IPD/XA-IR-00/3D		
			2	IPD/XA-IR-02/3D		
	RP	Engaging	3	IPD/XA-IR-04/3D	IPD/KT-IR-00	
			1	IPD/XA-IR-01/3D		
		Non-engaging	2	IPD/XA-IR-03/3D		IPD/KT-IR-01
			3	IPD/XA-IR-05/3D		
Medentis® ICX	RP	Engaging	1	IPD/YB-IR-00/3D	IPD/KT-IR-00	
			2	IPD/YB-IR-02/3D		
			3	IPD/YB-IR-04/3D		
		Non-engaging	1	IPD/YB-IR-01/3D		IPD/KT-IR-01
			2	IPD/YB-IR-03/3D		
			3	IPD/YB-IR-05/3D		